

New VariFit™ NIV Mask

Complementing our
existing range of critical
care patient interfaces



To arrange
a demonstration
call 0118 9656 376
www.intersurgical.co.uk/contact



Quality, innovation and choice

Interact with us



www.intersurgical.co.uk

VariFit™ NIV Mask

Introducing the VariFit NIV mask offering a further choice to meet your non-invasive ventilation needs. VariFit is a single patient use hospital NIV mask for adult patients featuring advanced AIRgel® cushion technology for enhanced patient comfort. The range includes three variations:

- Vented Mask with an anti-asphyxiation valve for use with single limb breathing systems that do not have a CO₂ leak port
- Non-vented Mask with an anti-asphyxiation valve for use with single limb breathing systems that have a CO₂ leak port
- Non-vented Mask for use with dual limb breathing systems

These variations can be identified by the colour of the mask elbow.

VariFit NIV masks are available in four sizes with a size guide to help select the appropriate mask.

For further information on our full range visit:
www.intersurgical.co.uk/products/critical-care/patient-interfaces

AIRgel is a registered trade mark of Sleepnet Corporation.
VariFit NIV masks are designed and manufactured by Sleepnet Corporation,
Hampton, NH, 03842 USA and distributed by Intersurgical Ltd, Crane
House, Molly Millars Lane, Wokingham, Berkshire, RG41 2RZ, UK



Quality, innovation and choice

T: +44 (0)118 9656 300
info@intersurgical.co.uk
www.intersurgical.co.uk



Also in our range of patient interfaces...



FaceFit™ NIV Mask



StarMed Respiratory Hoods



INTERSURGICAL
COMPLETE RESPIRATORY SYSTEMS



Please
recycle