

# Intersurgical's Universal Neonatal Resuscitation Breathing System

For use in both Delivery and Neonatal Units - Improving inventory management, reducing unnecessary waste.

10mm Flextube neonatal resuscitation breathing system with variable PEEP, double swivel elbow, Neopuff® and universal connectors, 1.2m




To view the neonatal resuscitation breathing systems with variable PEEP range, scan the QR code.

## 0.8m, 1.2m and 2m tubing options

Offering more choice for different ventilatory systems

## Click adjustment on PEEP cap

No accidental adjustment during operation



**Versatile connection**  
22M/15F fits all patient interfaces

Code	Description	Box Qty.
6432001	10mm Flextube neonatal resuscitation breathing system with variable PEEP, double swivel elbow for Neopuff®, 2m	10
6431000	10mm Flextube neonatal resuscitation breathing system with variable PEEP, double swivel elbow, 1.2m	15
6430000	10mm Flextube neonatal resuscitation breathing system with variable PEEP, double swivel elbow, 0.8m	15
 6433000	10mm Flextube neonatal resuscitation breathing system with variable PEEP, double swivel elbow, Neopuff® and universal connectors, 1.2m	15
6432000	10mm Flextube neonatal resuscitation breathing system with variable PEEP, double swivel elbow for Neopuff®, 1.2m	15
1661000	Straight connector 15M/11.5F - 10M	25



Quality, innovation and choice

If you want further information/support, or you wish to place an order, please contact me:  
John McRink - 07765643636 [jm@intersurgical.co.uk](mailto:jm@intersurgical.co.uk)



Interact with us  
[www.intersurgical.co.uk](http://www.intersurgical.co.uk)