

# 7 steps to improve patient and environmental safety within anaesthesia

## Integral monitoring line

Intersurgical's modular breathing systems include an integral gas monitoring line which helps to **improve accuracy** and **reduce patient connections**.



1

## Silver Knight™

Antimicrobial silver ions within our Silver Knight™ breathing systems help to **reduce bacteria by 99.9%** and help **minimise cross-contamination**.



2

## Tethered caps

All Intersurgical anaesthetic breathing systems have brightly coloured and tethered luer caps, which help **reduce the risk of foreign body aspiration during ventilation**.



3

## Push and Twist

Intersurgical's unique Push and Twist feature ensures a **secure, leak free fit every time**.



4

## Colour coded filters


Intersurgical's filter range is colour coded for **ease of identification** and to ensure the correct filter/HME is used, **minimising any potential harm to the patient**.



5

## Intersurgical sodalime

Fewer kilos of Spherasorb™ are used per hour of anaesthesia than many other absorbents and production waste is recycled back into usable soda lime, **reducing the amount of waste**.



6

## Reduce your environmental impact

Intersurgical's ClearLite™ anaesthetic mask offers a **72% reduction in environmental impact** compared to a traditional PVC mask.



7

