




# Intersurgical Venturi valves – Minimum flow rates




**Product Code: 1400024**  
Venturi valve 24% oxygen  
Blue  
2 L/min  
NPC Code: FDD718




**Product Code: 1400028**  
Venturi valve 28% oxygen  
White  
4 L/min  
NPC Code: FDD720



**Product Code: 1400035**  
Venturi valve 35% oxygen  
Yellow  
8 L/min  
NPC Code: FDD724



**Product Code: 1400040**  
Venturi valve 40% oxygen  
Red  
10 L/min  
NPC Code: FDD726



**Product Code: 1400060**  
Venturi valve 60% oxygen  
Green  
15 L/min  
NPC Code: FDD728