



Neonatal Resuscitation

with variable PEEP



Critical Care ▪ Breathing Systems for Critical Care

Neonatal resuscitation breathing systems with variable PEEP

A range of 10 mm Neonatal resuscitation products which incorporate a click adjustable orifice control on the patient elbow to provide variable PEEP.

The breathing systems are available in a variety of lengths to accommodate various types of equipment and positioning of the patient.

These breathing systems can be used in NICU and obstetric units with the following resuscitation platforms and gas drivers:

- Dräger®
- Viamed Tom Thumb™
- Fisher & Paykel® Neopuff®

A selection of tubing options

Offering more choice for different ventilatory systems

Click adjustment on PEEP cap


No accidental adjustment during operation



22M/15F swivel patient connection

Optimise tube positioning and patient comfort

Make an enquiry



Code	Description	Box Qty.
6430000	10 mm Flextube™ neonatal resuscitation breathing system with variable PEEP, double swivel elbow and 15F/10F/15M/11.5F adaptors. ≥ 0.8 m	15
6431000	10 mm Flextube neonatal resuscitation breathing system with variable PEEP, double swivel elbow and 15F/10F adaptors. ≥ 1.2 m	15
6432000	10 mm Flextube neonatal resuscitation breathing system for use with Neopuff® resuscitators with variable PEEP, double swivel elbow. ≥ 1.2 m	15
6431002	10 mm Flextube neonatal resuscitation breathing system for use with Neopuff® and universal resuscitators with variable PEEP, double swivel elbow. ≥ 2.0 m	15
6433000	10 mm Flextube neonatal resuscitation breathing system for use with Neopuff® and universal resuscitators with variable PEEP, double swivel elbow. ≥ 1.2 m	15

Illustrated

Dräger is a registered trademark of Drägerwerk AG & Co. KGaA. Tom Thumb is a trademark of Viamed, Ltd. Fisher & Paykel is a registered trademark of Fisher & Paykel Healthcare Limited. Neopuff is a registered trademark of Fisher & Paykel Healthcare Limited.

IS1.9 • Issue 10 10.24



Intersurgical Ltd, Crane House, Molly Millars Lane, Wokingham, Berkshire, RG41 2RZ, UK
T: +44 (0)118 965 6300 info@intersurgical.com www.intersurgical.com



The manufacturer Intersurgical Ltd is certified to ISO 14001:2015, ISO 9001:2015, ISO 13485:2016 and MDSAP

Please think before you print
Save energy and paper.
If you must print this information sheet please print it double sided.