



# Solus™ laryngeal mask airways

A comprehensive range of single use laryngeal mask airways



Airway Management • Airway Devices

## Solus™ laryngeal mask airways

A comprehensive range of single use laryngeal mask airways suitable for anaesthesia and emergency medicine.

Designed to meet the needs of a variety of clinical scenarios, all Solus laryngeal mask airways are single use, latex free and supplied sterile.

**Low friction material**  
Atraumatic insertion & removal

**Classic cuff shape**  
Renowned for reliability

**Essential information**  
Clearly visible for quick reference

**Integral inflation line**  
Reduces the risk of damage whilst in use

**High quality valve**  
To ensure continual cuff integrity

[▶ Watch the video](#)

### Solus Standard

Using specially chosen materials, the Solus Standard is manufactured with the quality, flexibility and strength required to work effectively in a range of demanding clinical applications, including paediatric anaesthesia. Available in a full size range, from adult to neonate.



[Make an enquiry](#)

Code	Description	Size	Weight	Box Qty.
8005000	Solus Standard, laryngeal mask airway	5 Large adult	70+kg	20
8004000	Solus Standard, laryngeal mask airway	4 Medium adult	50–70kg	20
8003000	Solus Standard, laryngeal mask airway	3 Small adult	30–50kg	20
8025000	Solus Standard, laryngeal mask airway	2.5 Large paediatric	20–30kg	20
8002000	Solus Standard, laryngeal mask airway	2 Small paediatric	10–20kg	20
8015000	Solus Standard, laryngeal mask airway	1.5 Infant	5–10kg	20
8001000	Solus Standard, laryngeal mask airway	1 Neonate	<5kg	20






## Solus Satin

Solus Satin is our range of single use laryngeal mask airways featuring a softer airway tube providing increased flexibility.



Make an enquiry

Code	Description	Size	Weight	Box Qty.
8005003	Solus Satin, laryngeal mask airway	5 Large adult	70+kg	20 
8004003	Solus Satin, laryngeal mask airway	4 Medium adult	50–70kg	20 
8003003	Solus Satin, laryngeal mask airway	3 Small adult	30–50kg	20 

## Solus MRI-safe

A range of single use, sterile, laryngeal mask airways fitted with specially designed non-ferrous valves, guaranteed not to interfere with the magnet in an MRI scanner. The plastic valve has been selected to ensure full reliability throughout the shelf life of each device.



Make an enquiry




Code	Description	Size	Weight	Box Qty.
8005002	Solus MRI-safe, laryngeal mask airway	5 Large adult	70+kg	20 
8004002	Solus MRI-safe, laryngeal mask airway	4 Medium adult	50–70kg	20 
8003002	Solus MRI-safe, laryngeal mask airway	3 Small adult	30–50kg	20 
8025002	Solus MRI-safe, laryngeal mask airway	2.5 Large paediatric	20–30kg	20 
8002002	Solus MRI-safe, laryngeal mask airway	2 Small paediatric	10–20kg	20 
8015002	Solus MRI-safe, laryngeal mask airway	1.5 Infant	5–10kg	20 
8001002	Solus MRI-safe, laryngeal mask airway	1 Neonate	<5kg	20 

### Solus Curve

Solus Curve is our laryngeal mask for those who prefer the insertion characteristics of a device with a curved airway tube.




Make an enquiry






Code	Description	Size	Weight	Box Qty.
8005050	Solus Curve, laryngeal mask airway	5 Large adult	70+kg	20 
8004050	Solus Curve, laryngeal mask airway	4 Medium adult	50–70kg	20 
8003050	Solus Curve, laryngeal mask airway	3 Small adult	30–50kg	20 

### Solus Flexible

Solus Flexible provides an ideal solution for shared airway procedures such as ENT, dental, oro-maxillary and eye surgery where the rigidity of airway devices can obscure surgical access. The wire-reinforced tube permits flexion without kinking and can be moved at any time mid-procedure without concern for the cessation of gas flow.



Make an enquiry

Code	Description	Size	Weight	Box Qty.
8005001	Solus Flexible, wire-reinforced laryngeal mask airway	5 Large adult	70+kg	20 
8004001	Solus Flexible, wire-reinforced laryngeal mask airway	4 Medium adult	50–70kg	20 
8003001	Solus Flexible, wire-reinforced laryngeal mask airway	3 Small adult	30–50kg	20 
8025001	Solus Flexible, wire-reinforced laryngeal mask airway	2.5 Large paediatric	20–30kg	20 
8002001	Solus Flexible, wire-reinforced laryngeal mask airway	2 Small paediatric	10–20kg	20 

 Sterile

IS6.1 ▪ Issue 15 06.24



Intersurgical Ltd, Crane House, Molly Millars Lane, Wokingham, Berkshire, RG41 2RZ, UK  
 T: +44 (0)118 965 6300 info@intersurgical.com www.intersurgical.com



The manufacturer Intersurgical Ltd is certified to ISO 14001:2015, ISO 9001:2015, ISO 13485:2016 and MDSAP

**Please think before you print**  
 Save energy and paper.  
 If you must print this information sheet please print it double sided.