



STARMed Ventukit

Hood for CPAP therapy with built-in Venturi flow driver



Critical Care ▪ Patient Interfaces ▪ **STARMed**

STARMed Ventukit CPAP hood with Venturi flow driver

StarMed Ventukit

StarMed Ventukit is a complete and easy to use system simply activated by connecting to the O₂ source available in the hospital departments. StarMed Ventukit unique design consists of oxygen tubes and a Venturi flow driver integral to a CPAP hood. The kit also includes an adjustable PEEP valve for immediate therapy activation.

Read more evidence on www.intersurgical.co.uk

Technical Features

- The integrated Venturi flow driver can generate high patient flows while consuming a reduced amount of compressed oxygen
- O₂% can be adjusted from 40% to 100% at each PEEP level
- The adjustable PEEP valve can provide patient pressures in a range of settings up to 20 cmH₂O
- The user is guided by a performance table integral to the product, for independent PEEP and O₂% settings

Ventukit CPAP hood with Venturi flow driver, bi-directional anti-asphyxiation valve and patient access port, adjustable PEEP valve, ear plugs and measuring tape

Performance table subject to changes. Always refer to product labelling.

Requirement	Setting	Flow delivered
O ₂ % approx	PEEP cmH ₂ O	Flow O ₂ L/min A + B
40%	5	9 + 0
	7.5	9 + 0
	10	10 + 0
	12.5	10 + 0
	15	11 + 0
50%	20	12 + 0
	5	8 + 9
	7.5	8 + 9
	10	9 + 9
	12.5	9 + 8
	15	10 + 8
	20	12 + 9

Requirement	Setting	Flow delivered
O ₂ % approx	PEEP cmH ₂ O	Flow O ₂ L/min A + B
60%	5	7 + 16
	7.5	8 + 16
	10	8 + 14
	12.5	9 + 13
	15	9 + 13
100%	20	10 + 14
	Feed with a minimum O ₂ flow of 40 L/min	

Code	Size	Neck circumference	Box Qty.
CV100XXL/2W	Extra Extra Large	≥50cm	5
CV100XL/2W	Extra Large	45-52cm	5
CV100L/2W	Large	40-47cm	5
CV100M/2W	Medium	34-41cm	5
CV100S/2W	Small	27-34cm	5
CV100XS/2W	Extra Small	17-27cm	5

www.intersurgical.co.uk/info/starmed
Make an enquiry

Dual flow meter manufactured by FLOW METER S.p.A.

Intersurgical Ltd, Crane House, Molly Millars Lane, Wokingham, Berkshire, RG41 2RZ, UK
T: 0118 965 6300 F: 0118 965 6356 info@intersurgical.co.uk www.intersurgical.co.uk

Management System
ISO 9001:2015
EN ISO 13485:2016

Manufacturer: Intersurgical S.p.A. Mirandola (MO) Italy
Please think before you print
 Save energy and paper. If you must print this information sheet please print it double sided.