



OptiLube™ sachets of lubricant

For use with supraglottic airways, endotracheal tubes and naso/oropharyngeal airways



Airway Management ▪ Airway Accessories



OptiLube™ sachets of lubricant

For use with supraglottic airways, endotracheal tubes and naso/oropharyngeal airways

To complement our wide range of airway devices, we are delighted to announce the addition of sachets of lubricant to our product portfolio as a service to our customers.

OptiLube™ sachets of lubricant contain 5g of water soluble lubricating jelly for use in the preparation of airway devices such as supraglottic airways, endotracheal tubes and naso/oropharyngeal airways. It can also be used to prime the gastric channel of second generation airway devices such as the i-gel® or lubricate suction and orogastric tubes.

OptiLube sachets are harmless to human tissue and will not injure rubber, metal appliances, instruments and natural or synthetic textiles.



Code	Description	Box Qty.
5469000	OptiLube™, sachets of lubricant, 5g	150

[Make an enquiry](#)

OptiLube is a trademark of Optimum Medical Solutions.

Note: Not all products are available in all markets, please contact your local office for details.

IS6.7 • Issue 5.1 08.26



Intersurgical Ltd, Crane House, Molly Millars Lane, Wokingham, Berkshire, RG41 2RZ, UK
T: +44 (0)118 965 6300 info@intersurgical.com www.intersurgical.com



The manufacturer Intersurgical Ltd is certified to
ISO 14001:2015, ISO 9001:2015, ISO 13485:2016 and MDSAP

UK • Australia • Benelux • Brazil • Canada • China • Colombia • Czech Republic • Denmark • France • Germany • India • Ireland • Italy
Japan • Lithuania • Malaysia • Mexico • Norway • Philippines • Portugal • Russia • South Africa • Spain • Sweden • Taiwan • Turkey • USA