

# nFlow™

## Designed to care

The smallest patients often need the most care. Safety and comfort are vital to help infants with their earliest steps in life.

### nFlow infant CPAP


The Intersurgical nFlow system offers a complete care solution for infant CPAP with established 'gas flip' technology, providing a low 'work of breathing' option as an alternative to ventilator CPAP.

nFlow is compatible with flow drivers designed for use with gas flip technology CPAP, providing slightly higher system pressure associated with this treatment.

For more information visit: [www.intersurgical.co.uk/info/nflow](http://www.intersurgical.co.uk/info/nflow)



**A**



**B**

### Bonnets



**C**

### Nasal masks



**D**

### Nasal prongs




## Ordering made easy

We've made ordering easy. To make a complete system you can choose one product from group A and B, and choose a patient interface from C and/or D. At any time you can order separate items according to your needs.

Bonnet size guide

Masks and prongs size guide



|          | Code                           | Description  |                       | Box Qty. |
|----------|--------------------------------|--|-----------------------|----------|
| <b>A</b> | <b>nFlow breathing systems</b> |  |                       |          |
|          | 4700000                        | nFlow infant CPAP heated wire breathing system, 1.6m                                 |                       | 12       |
|          | 4700001                        | nFlow infant CPAP heated wire breathing system, 1.8m                                 |                       | 10       |
|          | 4700016                        | nFlow infant CPAP heated wire breathing system with detachable monitoring line, 1.6m |                       | 10       |
|          | 4708000                        | nFlow silencer kit   |                       | 50       |
|          | 4700011                        | nFlow replacement bonnet ties  |                       | 60       |
| <b>B</b> | <b>nFlow bonnets</b>           |  | <b>Size</b>           |          |
|          | 4707009                        | nFlow, infant bonnet, navy   | size 9 (40 - 42 cm)   | 10       |
|          | 4707008                        | nFlow, infant bonnet, turquoise  | size 8 (38 - 40 cm)   | 10       |
|          | 4707007                        | nFlow, infant bonnet, orange   | size 7 (36 - 38 cm)   | 10       |
|          | 4707006                        | nFlow, infant bonnet, red  | size 6 (34 - 36 cm)   | 10       |
|          | 4707005                        | nFlow, infant bonnet, green  | size 5 (32 - 34 cm)   | 20       |
|          | 4707004                        | nFlow, infant bonnet, light orange   | size 4 (30 - 32 cm)   | 20       |
|          | 4707003                        | nFlow, infant bonnet, blue   | size 3 (28 - 30 cm)   | 20       |
|          | 4707002                        | nFlow, infant bonnet, yellow   | size 2 (26 - 28 cm)   | 20       |
|          | 4707001                        | nFlow, infant bonnet, dark grey  | size 1 (24 - 26 cm)   | 20       |
|          | 4707110                        | nFlow, infant bonnet, pink   | size 0 (22 - 24 cm)   | 20       |
|          | 4707100                        | nFlow, infant bonnet, grey   | size 00 (20 - 22 cm)  | 20       |
|          | 4707000                        | nFlow, infant bonnet, white  | size 000 (18 - 20 cm) | 20       |
| <b>C</b> | <b>nFlow nasal masks</b>       |  | <b>Size</b>           |          |
|          | 4716000                        | nFlow nasal mask   | extra large           | 30       |
|          | 4706000                        | nFlow nasal mask   | large                 | 30       |
|          | 4705000                        | nFlow nasal mask   | medium                | 30       |
|          | 4704000                        | nFlow nasal mask   | small                 | 30       |
| <b>D</b> | <b>nFlow nasal prongs</b>      |  | <b>Size</b>           |          |
|          | 4703000                        | nFlow nasal prong  | large                 | 30       |
|          | 4702000                        | nFlow nasal prong  | medium                | 30       |
|          | 4701000                        | nFlow nasal prong  | small                 | 30       |
|          | 4711000                        | nFlow nasal prong  | extra small           | 30       |